



# Metrosep A Supp 10 Guard/2.0

6.1020.600

The Metrosep A Supp 10 Guard/2.0 column reliably protects the A Supp 10 analytical microbore separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column without tools

## Phạm vi giao hàng 6.1020.600

Qt.	Order no.	Description
1 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		

