



# Metrosep C Supp 1 - 150/4.0

6.1052.420

The Metrosep C Supp 1 - 150/4.0 separation column is the column of choice for the determination of low concentrations of standard cations. Detection limits below a µg/L range are achieved due to low baseline noise after sequential suppression.

## Phụ kiện tùy chọn

| Order no.  | Description                 |
|------------|-----------------------------|
| 6.1052.500 | Metrosep C Supp 1 Guard/4.0 |

The Metrosep C Supp 1 Guard/4.0 contains the C-Supp-1 column material and is used to protect Metrosep C-Supp-1 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C Supp 1 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.

