



## Sample rack, 12 x 250 mL, (PVC)

6.2041.310

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 12 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.

Plastic: Polyvinyl chloride (PVC)

### Phụ kiện tùy chọn

Order no.	Description
6.1432.320	Sample beaker, clear glass, 250 mL, 1 piece
	Sample beaker, clear glass, 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020 and 6.2068.030, 1 piece



6.1432.323 Sample beaker, amber glass, 250 mL, 1 piece

Sample beaker, amber glass, 250 mL, suitable for sample racks:  
6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020, and 6.2068.030, 1 piece



---

6.1453.220 Sample beaker, plastic (PP), 200 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample  
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020, and 6.2068.030, 1 piece



---

6.1453.250 Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample  
racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,  
6.2068.020, and 6.2068.030, 1 piece



---

6.1458.010 Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.

With 6xSGJ 14 and 3x SGJ 9 openings.

