



Low Volume dialysis cell

6.2729.200

Including four 6.2744.060 stoppers. For rapid Inline Dialysis.

Phạm vi giao hàng 6.2729.200

Qt.	Order no.	Description
4 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Phụ kiện tùy chọn

Order no.	Description	
6.2714.010	Dialysis membrane (cellulose acetate) Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis	
6.2714.030	Dialysis membrane (polyamide) For inline dialysis and filtration	
6.2744.000	Pressure screw PVDF With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)	