



# OMNIS Sample Robot S – NIR

2.1073.0010


OMNIS Sample Robot S – NIR 8 – 28 mm 77

OMNIS Sample Robot S – NIR

## Scope of delivery 2.1073.0010

| Qt.   | Order no.   | Description                |
|-------|-------------|----------------------------|
| 1 PCS | 1.1073.0010 | OMNIS Sample Robot S – NIR |





OMNIS Sample Robot S – NIR for the direct transition to fully automated spectroscopy. The system provides space for 1 sample rack for up to 77 vials with diameters of 8–28 mm.  
The OMNIS Sample Robot S – NIR supports the analyses of both liquid and solid samples.



Function license for utilization of the command parameters for the  
OMNIS Sample Robot S – NIR.



## Optional accessories

| Order no.   | Description  |   |
|-------------|--|---|
| 6.02041.080 | <b>OMNIS sample rack 77 x 8 mm</b><br><br>OMNIS sample rack for OMNIS Sample Robot – NIR suitable for 77 disposable vials with an optical path length of 8 mm for measurements in transmission (6.7402.240). |    |
| 6.02041.090 | <b>OMNIS sample rack 77 x 15 mm</b><br><br>OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 77 disposable vials with a diameter of 15 mm for measurements in reflection (6.7402.110).            |    |
| 6.02041.100 | <b>OMNIS sample rack 50 x 19 mm</b><br><br>OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 50 disposable vials with a diameter of 19 mm for measurements in reflection (6.7402.120).            |  |
| 6.02041.110 | <b>OMNIS sample rack 50 x 22 mm</b><br><br>OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 50 disposable vials with a diameter of 22 mm for measurements in reflection (6.7402.130).            |  |

**6.02041.120**

**OMNIS sample rack 27 x 28 mm**

OMNIS sample rack for OMNIS Sample Robot – NIR, suitable for 27 disposable vials with a diameter of 28 mm for measurements in reflection (6.7402.140).

---

