



728

2.728.0040

645685665/765715/775725/776 DosimatTitrino6.2001.010

Scope of delivery 2.728.0040

Qt.	Order no.	Description
1 PCS	1.728.0010	728 Stirrer without stand
For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.		



1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2001.010

Support stand with support rod

For Dosimats and Titrinos



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2723.130

Cover for 655/665/765

Dust cover for 655, 665, 765 Dosimats

