



804

2.804.0010

802 802

Scope of delivery 2.804.0010

Qt.	Order no.	Description
1 PCS	1.804.0010	804 Ti Stand without stand rod
Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.		



Optional accessories

Order no.	Description	
6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrande	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug	