



# 826 pH mobile Primatrode

2.826.0020

pH IrDA NTC Primatrode 6.0228.020IP67

## Scope of delivery 2.826.0020

Qt.	Order no.	Description
1 PCS	1.826.0010	826 pH mobile

Portable pH meter with IrDa interface, including carrying strap and batteries.



1 PCS

**6.2050.000**

**Carrying strap for pH Meter**

Carrying strap for 704, 826 pH meter.



1 PCS

**6.2133.000**

**Batteries LR6, 1.5V, 4 pieces**

1.5 V LR6 type. Set of 4 pieces. Used with 704 pH Meter, 773 pH/mV Simulator




1 PCS

**6.2621.140**

**Hexagon key 2.5 mm**



## Optional accessories

Order no.	Description	
6.2307.230	<b>Buffer solutions pH 4, 7 and 9</b>  Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL	
6.2325.000	<b>pHit kit</b>  Maintenance kit for pH electrodes  The kit contains: <ul style="list-style-type: none"><li>• 50 mL cleaning solution</li><li>• 50 mL 3M KCl solution</li><li>• 50 mL storage solution</li><li>• 2 Storage vessels</li><li>• Instructions for use</li></ul>	