

pH Lab Irda 120V US

2.827.0015

Scope of delivery 2.827.0015

Qt.	Order no.	Description
1 PCS	6.2013.010	Clamping ring For support rods with a diameter of 10 mm.
		
1 PCS	6.2016.050	Support rod / 300 mm
		
1 PCS	6.2021.020	Electrode holder Electrode holder for 4 electrodes and 2 buret tips
		

1 PCS

6.2161.020

Plug-in power supply unit US 120V/6V DC

Power supply for 827 pH lab



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

