




IC Conductivity Detector MB

2.850.9020

DSP

-
- (2mm) IC-MS IC-ICP/MS
- 0 ... 15000 $\mu\text{S}/\text{cm}$
- 0.3 μL
- X2CrNiMo17-12-2 (316 L) MSA
- 10.0 MPa (100 bar)
- 20 ... 50 $^{\circ}\text{C}$ 5 $^{\circ}\text{C}$
- < 0.001 $^{\circ}\text{C}$
- < 0.2 $^{\circ}\text{nS}/\text{cm}$
- ID 0.18 mm

MagIC Net 4.1

Scope of delivery 2.850.9020

Qt.	Order no.	Description
1 PCS	6.2744.014	Pressure screw 2x With UNF 10/32 connection. For the connection of PEEK capillaries

