





900 Touch Control

2.900.0010

Titrand USB 856 Conductivity Module 867 pH Module 846 Dosing Interface USB U USB

Scope of delivery 2.900.0010

Qt.	Order no.	Description
1 PCS	1.900.0010	900 Touch Control
Operating unit for the Titrandos, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to Favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick. Dialog languages: German, English, Chinese, French, Spanish, Portuguese, Russian, Korean, Polish and Italian.		



1 PCS

6.2723.300

Protective cover for 900 Touch Control

Protective cover for 900 Touch Control



1 PCS

6.6064.010

USB Memory Stick for 900 Touch Control

