




Booster20A

3500003600


Booster20A PGSTAT302N 20 Booster20A 20 V

FRA32M Autolab Booster20A booster Booster20A DCAC

Scope of delivery 3500003600

Qt.	Order no.	Description
1 PCS	BST.EMER. BUTN	Booster emergency stop button

The Booster20A is fitted with an emergency stop button which can be position in close proximity to the instrument. This button can be used to immediatly disconnect the Booster20A from the cell if necessary.



1 PCS

CABLE.PWR Power cable

Standard power cable for Autolab instruments and accessories.



1 PCS

**CBL.BSTR20A. Booster20A digital interface cable
DIO**

Digital control cable for the Autolab Booster20A.

