



Pt1000 Ecotrode Gel

6.00221.600

pH Pt1000 pH cKCl= 6.2308.000

Scope of delivery 6.00221.600


Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.2104.600	Electrode cable for plug in head U/plug F, 2x2 mm B, 1m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2104.610	Electrode cable for plug in head U/plug F, 2x2 mm B, 2m For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL Electrolyte KCl (sat.), 250 mL	