



Metrosep A Supp 19 - 100/4.0

6.01034.410

- Metrosep A Supp 19

Metrosep A Supp 19 100 mm - 7

Metrosep A Supp 19 - 100/4.0

Scope of delivery 6.01034.410

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description
6.01034.500	Metrosep A Supp 19 Guard/4.0 <p>The Metrosep A Supp 19 Guard/4.0 reliably protects the anion separation columns of the Metrosep A Supp 19 product family against contaminations from the sample or eluent and thus prolong their service life significantly. The guard columns and separation columns of the Metrosep A Supp 19 product family are made of PEEK and filled with the same material. This ensures that the chromatographic separating efficiency is not restricted in any way.</p> <p>The "On Column Guard System" makes it possible to screw the guard column onto the anion separation column directly and conveniently. The economical price and simple handling make the Metrosep A Supp 19 Guard/4.0 the ideal protection for the separation column.</p>

