



# Metrosep A Supp 21 - 250/4.0

6.01036.430

Metrosep A Supp 21 /4.6 µm pH

Metrosep A Supp 21 - 250/4.0 DCAA US EPA 300.1 A DIN EN ISO 10304-1&4 µg/L

## Scope of delivery 6.01036.430

| Qt.   | Order no.  | Description      |
|-------|------------|------------------|
| 2 PCS | 6.2744.060 | Threaded stopper |

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.

