



## Metrosep C 6 - 250/2.0

6.01051.230

Metrosep C 6 - 250/2.0 C-6

IC-MS

### Scope of delivery 6.01051.230

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Optional accessories

Order no.	Description
6.01051.600	<b>Metrosep C 6 Guard/2.0</b> <p>The Metrosep C 6 Guard/2.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/2.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.</p>

