



3 NS14/ 4 6.4mmP&P

6.01403.030

3 NS14 4 6.4 mm OMNIS Sample Robot Pick&Place 2 Pick&Place 120 mL

Scope of delivery 6.01403.030

Qt.	Order no.	Description
2 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



2 PCS

6.1446.000

SGJ stopper / B-14/(15)



2 PCS

6.1808.000

Connection bushing M6



2 PCS


6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Optional accessories

Order no.	Description	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	 A photograph showing three replacement tubings. Two are light grey and one is dark red. They are arranged in a slightly overlapping, curved manner, showing their tapered ends.
