



OMNIS Polytron 601406XXX

6.01405.040

Polytron OMNIS Titrator ·

Scope of delivery 6.01405.040

Qt.	Order no.	Description
1 PCS	6.1454.130	Form seal made of FPM for 61414xxx for KF and VA instruments



Optional accessories

Order no.	Description
-----------	-------------

6.02709.020 **Septum stopper**

Plastic septum stopper, suitable for 6014050X0 titration vessel lids.



6.02709.030 **Screw nipple**

Plastic screw nipple, suitable for 6014050X0 titration vessel lids.



6.02709.040 **Stopper**

Stopper for titration vessel lids 6014050X0.

