



20 – 90 mL

6.01406.220

6014050X0 20 – 90 mL

Optional accessories

Order no.	Description
6.01405.000	5 x SGJ 14 titration vessel lid OMNIS

Titration vessel lid with five SGJ 14 openings for 601406XX0 titration vessels with screw cap.



6.01405.010

OMNIS titration vessel lid for volumetric Karl Fischer titration

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator.



6.01405.040

OMNIS titration vessel lid for 601406XXX titration vessels for Polytron

Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator with Polytron.

