




5072 mm

6.02601.020

OMNIS Sample Robot Pick&Place 5072 mm

Optional accessories

Order no.	Description
6.01400.000	Sample beakers, clear glass, 250 mL, 10 pieces Sample beakers, clear glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.



6.01400.003 Sample beakers, amber glass, 250 mL, 10 pieces

Sample beakers, amber glass, 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.



6.01400.100 Sample beakers, plastic (PP), 250 mL, 10 pieces

Sample beakers, plastic (polypropylene, PP), 250 mL, suitable for 602041010 OMNIS sample rack, 10 pieces.



6.02064.020 Retainer for 200 mL sample beaker

Retainer for fixing the 200 mL sample beaker in place in the gripper fingers of the OMNIS Sample Robot Pick&Place.



6.1459.310 Sample beaker, disposable, plastic (PP), 200 mL, 1000 pieces

Sample beaker, disposable, plastic (polypropylene, PP), 200 mL, suitable for 6.02041.020 OMNIS Software sample rack, 6.2041.370 und 6.2041.830 sample racks, 1000 pieces.

