




KF Dis-Cover 120 mL OMNIS 16

6.02710.010

OMNIS Sample Robot Pick&Place 6.01400.2006.01400.300 6.01400.303 ·

Scope of delivery 6.02710.010

Qt.	Order no.	Description
1 PCS	6.05700.170	Consumable Kit seals for the 6.02710.010 KF Dis-Cover lid

Replacement seals for 6.02710.010 lid, 16 pieces



Optional accessories

Order no.	Description	
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.300	Sample beakers, clear glass, 120 mL, 20 pieces Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.303	Sample beakers, amber glass, 120 mL, 20 pieces Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces	