



2

6.02729.110

858 Professional Sample Processor / 919 IC Autosampler plus

Scope of delivery 6.02729.110

Qt.	Order no.	Description
1 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
6.02057.030	Ultrafiltration cell 2 holder Holder for mounting the 6.02729.110 ultrafiltration cell 2 on the 858 Professional Sample Processor / 919 IC Autosampler plus.	
6.02057.120	Ultrafiltration cell 2 holder, stand-alone Holder for the 6.02729.110 ultrafiltration cell 2	
6.2714.020	Filtration membrane reg. cellulose, 0.2 µm, 50 pieces Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.	
6.2714.030	Dialysis membrane (polyamide) For inline dialysis and filtration	

6.2714.050

Filtration membrane reg. cellulose, 0.45 µm, 50 pieces

Nominal pore size: 0.45 µm. Set of 50 pieces. For inline filtration.



6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



6.2831.010

Plastic tweezers

For dialysis and filtration membranes

