



6.0301.100

### Scope of delivery 6.0301.100


Qt.	Order no.	Description
1 PCS	6.1236.050	<b>Sleeve with SGJ 14/12 mm</b> Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



# Optional accessories

Order no.	Description	
6.2106.020	<b>Strand / 1 m / 2 x B</b>  Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2106.060	<b>Strand / 2 m / 2 x B</b>  Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	