



iAg

6.0450.300

pH

-
-
-
-
-

$c(\text{KNO}_3) = 1 \text{ mol/L}$

iTrodes TitrandTi-Touch 913/914 Metern

Scope of delivery 6.0450.300

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description
-----------	-------------

2.854.0010 854 iConnect

The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).



6.2310.010 Electrolyte KNO3-1mol/L 250 mL

Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)

