



Ca

6.0510.100

ISE

- Ca^{2+} $5 \cdot 10^{-7}$ 1 mol/L
-

/

$\text{c}(\text{NH}_4\text{NO}_3) = 1 \text{ mol/L}$

Scope of delivery 6.0510.100

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2327.000	Electrolyte NH4NO3 1 mol/L (50 mL) Electrolyte solution $c(\text{NH}_4\text{NO}_3) = 1 \text{ mol/L}$ (for combined Ca-ISE 6.0510.100)	