



6.05800.030

2 OMNIS Sample Robot Pick&Place 4

Scope of delivery 6.05800.030

Qt.	Order no.	Description
1 PCS	6.01543.100	Aspiration tip automation

Aspiration tip for automated OMNIS titrations.




3 PCS

6.01805.000

FEP tubing connection / M6 / 66.5 cm

With light protection and protection against kinks.



1 PCS

6.01805.010

PTFE tubing connection / M8 / 66.5 cm

With light protection and protection against kinks.



1 PCS

6.05001.010

**"Peristaltic" (2-channel) pump module
accessory set**

The accessory set contains all of the components required for the tube connections between a rinsing pump and an aspiration pump at the "Peristaltic" pump module. In addition to preassembled tubing with Luer connectors, each of these also includes two press clamps and pump tubing.



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1621.100

Cap for 10-L container

Cap for 10-L container. For easy and safe transportation of filled containers.



2 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.

