



Metrosep A Supp 4 - 250/4.0

6.1006.430

Metrosep A Supp 4 - 250/4.0 9 µm Metrosep A Supp 4 - 250/4.0 12.5 Metrosep Anion Dual 2 - 75/4.6A Supp 4 - 250/4.0

Optional accessories

Order no.	Description
6.01021.500	Metrosep A Supp 4 Guard/4.0

The Metrosep A Supp 4 Guard/4.0 reliably protects the Metrosep A Supp 4 - 250/4.0 anion columns against contamination from the sample or eluent.



The Metrosep A Supp 4 S-Guard/4.0 reliably protects the Metrosep A Supp 4 IC anion columns against contamination from the sample or eluent.

