



Metrosep C 3 - 250/4.0

6.1010.430

« » Metrosep C 3 - 250/4.0 42000 51000 31000

Metrosep C 3 - 250/4.0 Metrosep C 3 - 250/4.0

Optional accessories

Order no.	Description
6.1010.450	Metrosep C 3 Guard/4.0 The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C 3 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.



6.1010.460

Metrosep C 3 S-Guard/4.0

The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column.

