



Metrosep A Supp 10 - 250/4.0

6.1020.030

Metrosep A Supp 10 - 250/4.0 / 4.6 μ mA-Supp-10 MSM-HC

Metrosep A Supp 10 - 250/4.0

Scope of delivery 6.1020.030

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
6.1020.500	Metrosep A Supp 10 Guard/4.0 The Metrosep A Supp 10 Guard/4.0 reliably protects the Metrosep A Supp 10 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column. No tools are required.	
6.1020.510	Metrosep A Supp 10 S-Guard/4.0 The Metrosep A Supp 10 S-Guard/4.0 reliably protects the Metrosep A Supp 10 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	
6.1020.520	Metrosep A Supp 10 Guard HC/4.0 The Metrosep A Supp 10 Guard HC/4.0 is the high-capacity version of the Metrosep A Supp 10 Guard/4.0.	