



Metrosep A Supp 10 - 75/2.0

6.1020.270

Metrosep A Supp 10 - 75/2.0 / 4.6 µm2 mm Metrosep A Supp 10 4 mm

Metrosep A Supp 10 - 75/2.0 IC-MS

Scope of delivery 6.1020.270

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description
6.1020.600	Metrosep A Supp 10 Guard/2.0 The Metrosep A Supp 10 Guard/2.0 column reliably protects the A Supp 10 analytical microbore separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column without tools

