



Metrosep A Supp 16 - 250/4.0

6.1031.430

Metrosep A Supp 16 Metrosep A Supp 16 / Metrosep A Supp 16

Metrosep A Supp 16 - 250/4.0 Metrosep A Supp 16 - 250/4.0

Scope of delivery 6.1031.430

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
6.1031.500	Metrosep A Supp 16 Guard/4.0 The Metrosep A Supp 16 Guard/4.0 reliably protects the Metrosep A Supp 16 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.	
6.1031.510	Metrosep A Supp 16 S-Guard/4.0 The Metrosep A Supp 16 S-Guard/4.0 reliably protects the Metrosep A Supp 16 analytical separation columns against contamination. The guard column is connected to the separation column simply using capillary connections. No tools are required.	