



CVS Pt1000 9 cm

6.1110.120

Pt1000 B Professional CVS

Scope of delivery 6.1110.120

Qt.	Order no.	Description
1 PCS	6.1236.020	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm and O-ring.



Optional accessories

Order no.	Description
6.2104.140	Electrode cable / 1 m / 2 x 2 mm For connecting temperature sensors with Metrohm plug-in head G with 780, 781, 826, 827 pH/Ion meters and all Titrandos 