



O2-Lumitrode

6.1116.000

(DO) 913 pH/DO Meter 914 pH/DO DO

-
-
-
-

O₂

Optional accessories

Order no.	Description
-----------	-------------

6.1260.000

Calibration vessel O2

Calibration vessel for calibrating the O₂-Lumitrode at 100% air saturation.



6.2329.000

Oxygen standard 0%

Oxygen standard for 0% calibration of the O₂-Lumitrode. Contains 5 sachets of 30 mL each.



6.5623.000

O2 Lumitrode replacement set

Replacement set for the O₂ Lumitrode. The set contains a certified O₂ cap as replacement, a calibration vessel, and 5 x 30 mL oxygen standard, 0%.

