



# scTRACE Gold

6.1241.080

scTRACE Gold scTRACE Gold 6.1258.000 Professional VA VA VA Computrace 6.2112.130  
 scTRACE Gold scTRACE Gold

## Scope of delivery 6.1241.080


Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



## Optional accessories

Order no.	Description	
6.1258.000	<b>scTRACE Gold (set with 4 pieces)</b>  Gold micro-wire electrode on thick film basis for the determination of arsenic in waters. The scTRACE Gold is a combined voltammetric sensor for voltammetry. Gold micro-wire working electrode, Ag/AgCl reference electrode and carbon auxiliary electrode are located on a sensor. Simple in its handling and practically maintenance-free, the scTRACE Gold can for example be used to determine arsenic in water in the trace range. Together with the 6.1241.080 electrode shaft, it can be used with any Metrohm voltammetry measuring stand.	
6.2112.130	<b>Adapter plug contact 5.5 mm / screw contact, gold-plated</b>  For VA electrodes	