



Polytron 786

6.1462.250

Polytron 786 Polytron

Scope of delivery 6.1462.250

Qt.	Order no.	Description
1 PCS	6.1808.000	Connection bushing M6



Optional accessories

Order no.	Description	
6.1829.050	Replacement tubings for retracted tubings in titration heads and robotic arms Three replacement tubings with tapering tip.	
6.9012.000	Dispersing aggregate for Polytron, 125 mm Standard dispersing aggregate for Polytron 1300 D	
6.9012.010	Dispersing aggregate for Polytron, 157 mm Special dispersing aggregate (Metrohm version) with protruding knives	