



807 Dosing Unit ETFE 5 mL

6.1575.150

Dosino (HF) 807 Dosing Unit ETFE 5 mL

Scope of delivery 6.1575.150

Qt.	Order no.	Description
1 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Optional accessories

Order no.	Description	
6.1566.150	Cylinder element ETFE 5 mL For dosing units used for dosing alkali reagents or hydrofluoric acid (HF) with Dosinos. Less precision and reproducibility is to be expected in comparison with glass cylinders.	
6.1574.150	Cylinder element glass 5 mL Glass cylinder element complete with cylinder, piston, cylinder base, valve disk, and centering tube. For 807 Dosing Unit 5 mL.	