



6.5337.200

## Scope of delivery 6.5337.200

Qt.	Order no.	Description
1 PCS	6.1257.820	50 µm spacer to Thin-Layer cell

50 µm spacer for the amperometric Thin-Layer cell, 3 pcs.



1 PCS

**6.1257.840**

**25 µm spacer to Thin-Layer cell**

25 µm spacer for the amperometric Thin-Layer cell, 3 pcs.



---

1 PCS

**6.1831.010**

**PEEK capillary 0.25 mm i.d. / 3 m**

For all IC components.



---

1 PCS

**6.2621.140**

**Hexagon key 2.5 mm**



---

1 PCS

**6.2744.014**

**Pressure screw 2x**

With UNF 10/32 connection. For the connection of PEEK capillaries



1 PCS

**6.2802.000**

**Polishing set for solid-state electrodes**

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)

