



NIRS 10 100

6.7402.030

- NIRS DS2500 Analyzer2.922.0010
- NIRS XDS MasterLab Analyzer2.921.1310
- NIRS XDS MultiVial Analyzer2.921.1120
- NIRS XDS RapidContent Analyzer2.921.1110
- NIRS XDS RapidContent Analyzer - Solids2.921.1210

Scope of delivery 6.7402.030

Qt.	Order no.	Description
-----	-----------	-------------

The disposable backs are used for the secure sealing of the mini sample cups. The sample cups are used in the following:

- NIRS DS2500 Analyzer (**order number: 2.922.0010**)
- NIRS XDS MasterLab Analyzer (**order number: 2.921.1310**)
- NIRS XDS RCA Analyzer (**order number: 2.921.1110**)
- NIRS XDS MultiVial Analyzer (**order number: 2.921.1210**)



The package is comprised of 100 disposable backs.
