



6.9011.020

859 Titrotherm
Thermoprobe
HF

- (TAN) ASTM D8045
- (TBN)
-
- Ca/Mg
-

Scope of delivery 6.9011.020

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.

