



Booster 10A

BSTR10A.S

Booster10A PGSTAT100N, PGSTAT128N, PGSTAT302N, PGSTAT204 M204 10
Booster10A 20 V

FRA32M Autolab Booster10A booster Booster10A DCAC

Scope of delivery BSTR10A.S

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS	CABLE.PWR	Power cable
-------	-----------	-------------

Standard power cable for Autolab instruments and accessories.



Optional accessories

Order no.	Description	
AUTSRV. BSTR10.CELL.S	Booster10A Dummy cell	
	100 mOhm dummy cell for testing the Booster10A.	
