




EQCM.S

EQCM EQCM

sub $\mu\text{g}/\text{cm}^2$ EQCM 6 MHz AT

EQCM 6 MHz

Scope of delivery EQCM.S

Qt.	Order no.	Description
2 PCS	EQCM.AU	EQCM 6 MHz Au/TiO ₂ crystal working electrode
		EQCM 6 MHz Au/TiO ₂ crystal working electrode



1 PCS

EQCM.CE

EQCM counter electrode

Gold coil counter electrode for Autolab EQCM cell.



1 PCS

EQCM.REF.EL.S **Reference Electrode Ag/AgCl 3M KCl**

Reference Electrode Ag/AgCl 3M KCl

